

SKINTACT Diagnostic Resting ECG Electrodes

"SKINTACT Easitabs are designed for Resting 12 lead diagnostic ECG's. Easitabs are disposable eliminating the spread of cross infection associated with older style suction cups. As the name suggests Easitabs are easy to use, provide reliable high quality recordings, and are suitable for all types of patients."



Manufactured using Leonhad Lang's unique C-LINE technology which ensures that every electrode will be fresh, providing you with exceptional recording consistency at very affordable prices.

SKINTACT®

www.skintact.com



SKINTACT Diagnostic Resting ECG Electrodes

The SKINTACT Easitab System –

Your best choice for an excellent and safe resting ECG:

Excellent Trace Quality. The conductive adhesive Aqua-Tac hydrogel of SKINTACT Easitabs provides a reliable signal tracing of consistent high quality. Easitab F-WA00 and FS-WB00 with Aqua-Wet gel complete the range for demanding conditions.

Secure Adhesion and Patient Comfort. SKINTACT Easitabs feel soft but stick well to the patient's skin. Their contact design prevents lift-off. After the ECG they can be removed comfortably – without loss of hair.

Superior Safety. Proper use of SKINTACT Easitabs eliminates cross infection whereas suction electrode bred bacteria and fungi.

Time Efficiency. Unlike suction electrodes, SKINTACT Easitabs on a handy card are applied quickly on all types of patients. Aqua-Tac gel requires no clean up after the ECG. Easiprep and Easiclip tab adaptors further facilitate your success.

As an alternative we offer you a choice of SKINTACT F-RG with standard stud connector and SKINTACT FS-RB4 with adaptor for a 4 mm banana connector.

	Order Code	Gel, Connector	Backing, Adhesive	Size, Dimensions	Packaging
	RT-14	Aqua-Tac, Easitab	Paper/PET, Gel	14 x 34 mm, pediatric, card of 10	100/pouch 500/box 5.000/case
	RT-34	Aqua-Tac, Easitab	Paper/PET, Gel	22 x 34 mm, adult, card of 10	100/pouch 500/box 5.000/case
	RT-74	Aqua-Tac, Easitab	Paper/PET, Gel	26 x 34 mm, adult, card of 10	100/pouch 500/box 5.000/case
	F-WA00	Aqua-Wet, Easitab	Foam, Standard	32 x 46 mm, adult, card of 10	50/pouch 500/box 3.000/case
	FS-WB00	Aqua-Wet, Easitab	Foam, Enhanced	28 x 37 mm, adult, card of 10	100/pouch 500/box 3.000/case
	F-RG	Aqua-Wet, Snap	Foam, Standard	32 x 41 mm, universal, card of 6	30/pouch 1.200/box 4.800/case
	FS-RB4	Aqua-Wet, 4 mm Adaptor	Foam, Enhanced	28 x 44 mm, adult, card of 3	30/pouch 900/box 3.600/case

				
Product	Easiclip	Easiprep	RAZOR	59316
Packaging Standard	10/pouch	500/pouch	50/box	10/pouch

SKINTACT Diagnostic Resting ECG Electrodes meet the ANSI/AAMI standard EC12 for Disposable ECG Electrodes.