

SKINTACT Low Cost Multipurpose ECG Electrodes

"SKINTACT Low Cost Multipurpose ECG Electrodes are designed to provide a high quality trace for all ECG applications. The latest in technological and production methodology enables us to convert these technical advances into a new generation of high quality, low cost, multifunctional ECG electrodes."



Manufactured using Leonhard Lang's unique C-LINE technology which ensures that every electrode will be fresh, providing you with exceptional recording consistency at very affordable prices.

SKINTACT®



www.skintact.com



SKINTACT Low Cost Multipurpose ECG Electrodes

Raising the Standards in ECG technology.

Our SKINTACT ECG electrodes on cards of six are our no compromise attitude to quality and variety at very competitive prices:

Quick'n'Easy. SKINTACT F-RG, F-RG1, F-TC1, F-TF, FS-RG, FS-RG1, FS-TC1 and FS-TF come on handy cards of six pieces. A perforation in the middle can be used to tear the card apart giving you two strips of three.

Superior Design. The label material overlaps the electrode and creates a lift tab which allows the clinical staff to handle the electrode easily even when wearing surgical gloves.

Secure Adhesion. The enhanced adhesion of SKINTACT FS-RG, FS-RG1, FS-TC1 and FS-TF assures that the electrodes stick very securely ensuring good contact even when the patient is perspiring.

SKINTACT F-RG, F-RG1, F-TC1 and F-TF are all produced with standard adhesive to the same high quality as the FS range. Suitable for times when the enhanced adhesive is not required, particularly where the patient's skin is more sensitive.



Order Code	F-RG	F-RG1	F-TC1	F-TF
Size	Universal	Universal	Universal	Adult
Gel, Connector	Aqua-Wet, Snap	Aqua-Tac, Snap	Aqua-Tac, Snap	Aqua-Wet, Snap
Backing, Adhesive	Foam, Standard	Foam, Standard	Foam, Standard	Foam, Standard
Dimensions	32 x 41 mm	32 x 41 mm	35 x 41 mm	43 x 51 mm
Suggested Use	Monitoring, Resting ECG	Monitoring, ICU	Monitoring, ICU	Monitoring, ICU, CCU
Packaging Standard	30/pouch 1.200/box 4.800/case	30/pouch 1.500/box 6.000/case	30/pouch 1.500/box 6.000/case	30/pouch 1.200/box 4.800/case

Alternative Alternative Alternative Alternative

Order Code	FS-RG	FS-RG1	FS-TC1	FS-TF
Size	Universal	Universal	Universal	Adult
Difference	Foam, Enhanced	Foam, Enhanced	Foam, Enhanced	Foam, Enhanced
Suggested Use	Monitoring, ICU, CCU	Monitoring, ICU	Monitoring, ICU	Monitoring, ICU, CCU, Stress Test

All SKINTACT Multipurpose ECG electrodes meet the ANSI/AAMI standard EC12 for disposable ECG electrodes.